


quantitative

- Continuous, quantitative HER2 score captures a wide range of expression
- Interpreting the HER2 Result:
 - Negative: < 10.7
 - Equivocal: ≥ 10.7 to < 11.5
 - Positive: ≥ 11.5

Single Gene scores



HER2

The Oncotype DX® Breast Cancer Assay:

- Predicts the magnitude of chemotherapy benefit
- Quantifies the likelihood of breast cancer recurrence
- Provides clinical experience for multiple patient populations

oncotype DX® Breast Cancer Assay

www.oncotypeDX.com
 (866) ONCOTYPE • (866) 662-6897
 customerservice@genomichealth.com

Genomic Health is a CLIA-certified, high-volume reference laboratory.



Genomic Health, Oncotype DX and Recurrence Score are registered trademarks of Genomic Health, Inc. ©2008 Genomic Health, Inc. All rights reserved. GHI284 09/08



Quantitative HER2



oncotype DX® Breast Cancer Assay

Accurate, Precise, Reproducible.



Concordance with IHC:

The Oncotype DX[®] Breast Cancer Assay Report to Include a Quantitative HER2 Score

Although the Recurrence Score remains the best tool for assessing prognosis and prediction of chemotherapy benefit, HER2 is an important marker for therapeutic decision-making for patients with breast cancer, and its measurement may significantly impact the chosen course of treatment.

Quantitative HER2 Information:

The HER2 Score by the Oncotype DX Breast Cancer Assay is another measure for clinicians and patients to determine HER2 status, and it provides further insight into individual patients' breast cancer tumor biology together with the Recurrence Score.

Highly Concordant with IHC and FISH:

The HER2 score's concordance with IHC and FISH meets or exceeds the ASCO/CAP guidelines requirement of 95% for determination of HER2 status.

Added Clinical Information:

Many clinicians who rely on the Oncotype DX Breast Cancer Assay for treatment planning have requested that Genomic Health also report the quantitative HER2 Score as this is already a component of the Recurrence Score.

– To provide further clarification especially in cases of equivocal IHC and/or FISH results or discordance between FISH and IHC.

Addressing a Limitation with Current Testing:

The impact of preanalytical variability can be minimized by "normalization" strategies used in quantitative gene expression assessment as performed by quantitative RT-PCR by Oncotype DX.


Oncotype DX Assay and central laboratory

IHC: **95%** (CI: 92%,96%)

HER 2 Central IHC* vs Oncotype DX
(Current ASCO/CAP Guidelines)

	IHC +	IHC -	Total
RT-PCR by Oncotype DX +	94 (78%)	4 (1%)	98
RT-PCR by Oncotype DX -	27 (22%)	439 (99%)	466
Total	121	443	564

Concordance: 95%
95% CI (92%,96%)



- 27 (22%) IHC 3+ cases were negative by RT-PCR
- 4 (1%) IHC negative cases were positive by RT-PCR

* IHC analysis using HercepTest[®]

- E2197 Study¹
- N = 755 patients
- The overall range of HER2 expression is approximately 1,000-fold
- Range of HER2 expression for HER2 positive is approximately 16-fold

HER2

¹ Baehner FL et al. HER2 Concordance Between Central Laboratory Immunohistochemistry & Quantitative Reverse Transcription Polymerase Chain Reaction in Intergroup Trial E2197. Presented at the ASCO Breast Cancer Symposium, September 6, 2008; Washington, D.C. Abstract #13

Concordance with FISH:


Oncotype DX Assay and central laboratory

FISH: **97%** (CI: 96%,99%)

HER 2 Central FISH* vs Oncotype DX
(Current ASCO/CAP Guidelines)

	FISH +	FISH -	Total
RT-PCR by Oncotype DX +	55 (98%)	11 (3%)	66
RT-PCR by Oncotype DX -	1 (2%)	408 (97%)	409
Total	56	419	475

Concordance: 97%
95% CI (96%,99%)



- 1 case (2%) was FISH positive but negative by RT-PCR
- 11 cases (3%) were RT-PCR positive but FISH negative

* FISH analysis using Vysis

- Kaiser Study²
- N = 568 patients
- The overall range of HER2 expression is approximately 500-fold
- Range of HER2 expression for HER2 positive is approximately 16-fold

HER2

² Baehner FL et al. HER2 Assessment in a Large Kaiser Permanente Case-Control Study: Comparison of FISH and Quantitative RT-PCR Performed by Central Laboratories. Presented at the ASCO Breast Cancer Symposium, September 6, 2008; Washington, D.C. Abstract #41